pt Mepilex® Border Ag

Curativo autoaderente, antimicrobiano, composto por espuma de silicone suave

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
www.molnlycke.com/virtual-patent-marking

en Mepilex® Border Ag

Self-adherent antimicrobial soft silicone foam dressing

℞ only **Caution: Federal [US] Law restricts this device to sale by or on the order of a physician [or properly licensed practitioner].**

Product description

Mepilex Border Ag consists of a Safetac® soft silicone wound contact layer, an absorbent polyurethane foam pad containing silver sulfate and colored by activated carbon, a layer with super absorbent polyacrylate fibres, a non woven and breathable and waterproof film.

Mepilex Border Ag is a soft silicone foam dressing that absorbs exudate, maintains a moist wound healing environment and has antimicrobial properties.

Mepilex Border Ag has been shown to inactivate wound related pathogens up to 7 days *in vitro*.

Safetac® is a proprietary adhesive technology developed by Mölnlycke that minimises pain to patients and trauma to intact skin or wounds.

Dressing material content:

Polyurethane, silicone, polyacrylate, cotton, viscose, polyester, polyolefin and silver sulfate in the wound pad area which equals to 1.2 mg/cm^2 silver

Indications for use

Mepilex Border Ag dressing is indicated for the management of exuding wounds such as leg and foot ulcers, pressure ulcers, traumatic and surgical wounds, superficial and partial thickness burns. Mepilex Border Ag can also be used under compression bandaging.

Silver sulfate present in the dressing helps reduce microbial colonization on the dressing.

Precautions

- Do not use Mepilex Border Ag on patient and/or user with known hypersensitivity to the ingoing materials/components of the product.
- Clinicians/ healthcare professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Frequent or prolonged use of this product may result in transient or permanent discoloration of skin.
- In the event of clinical infection Mepilex Border Ag does not replace the need for systemic therapy or other adequate infection treatment.
- Prior to commencing radiation therapy remove Mepilex Border Ag if product is present in the treatment area. A new dressing can be applied following treatment.
- Avoid contact with electrodes or conductive gels during electronic measurements, e.g. electrocardiograms (ECG) and electroencephalograms (EEG).
- Do not use Mepilex Border Ag together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- Other than saline solution or water, the interaction of cleansing agents in combination with Mepilex Border Ag has not been demonstrated.
- The interaction of Mepilex Border Ag with other topical treatments has not been demonstrated.
- For external use only.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.
- Do not use after expiry date. If the product is used after the expiry date product properties cannot be ensured.

Instructions for use

Mepilex Border Ag should be used by or under the supervision of a qualified health care professionals.

Note that local hygiene procedures should be followed prior to and following the dressing change.

1. Cleanse the wound with saline solution or water according to standard clinical practice.
2. Dry the surrounding skin thoroughly.
3. Select appropriate dressing size. For best clinical outcomes, ensure that the absorbent pad covers the wound bed and overlaps the wound edges.
4. Remove the release films and apply the adherent side to the wound. Do not stretch

Each Mepilex Border Ag dressing may be left in place for up to 7 days, depending on the patient, the condition of the wound and surrounding skin, or as indicated by accepted clinical practice.

Mepilex Border Ag is intended for short-term use up to 4 weeks. For long-term use, a clinical assessment by a physician is recommended.

Special storage conditions and handling conditions

Mepilex Border Ag should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. Note that any color variances in Mepilex Border Ag do not affect product effectiveness or safety.

Disposal should be handled according to local environmental procedures.

Other information

If any serious incident has occurred in relation to the use of Mepilex Border Ag, you should report to Mölnlycke Health Care and to your local competent authority.

Mepilex and Safetac are registered trademarks of Mölnlycke Health Care AB.

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Mepilex y Safetac son marcas registradas de Mölnlycke Health Care AB.