pt Mepilex® Ag

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


Mepilex® Ag

Antimicrobial soft silicone foam dressing

℞ only **Caution:** Federal (US) Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Product description

Mepilex Ag consists of:

- a Safetac® soft silicone wound contact layer layer
- a flexible absorbent pad of grey polyurethane foam containing silver sulfate and colored by activated carbon
- a breathable and waterproof outer film

Mepilex Ag contains silver sulphate that releases silver ions to create an effective bacterial barrier and inactivates a wide range of wound related pathogens (bacteria and fungi), shown *in vitro*. By reducing the number of microorganisms, Mepilex Ag may also reduce odour.

Mepilex Ag has also been shown to inactivate wound related pathogens, up to 7 days *in vitro*.

Safetac® is a proprietary adhesive technology developed by Mölnlycke that minimises pain to patients and trauma to intact skin or wounds.

Dressing material content:

Polyurethane, silicone and silver sulfate which equals to 1.2 mg/cm^2 silver

Indications for use

Mepilex Ag dressing is indicated for the management of low to moderately exuding wounds such as leg and foot ulcers, pressure ulcers and partial thickness burns.
Silver sulfate present in the dressing helps reduce microbial colonization on the dressing.

Precautions

- Do not use Mepilex Ag on patient and/or user with known hypersensitivity to the ingoing materials/components of the product.
- Clinicians/healthcare professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.
- Prior to commencing radiation therapy remove Mepilex Ag if product is present in the treatment area. A new dressing can be applied following treatment.
- Do not use Mepilex Ag together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- For external use only.
- Frequent or prolonged use of this product may result in transient or permanent discoloration of skin.
- In the event of clinical infection Mepilex Ag does not replace the need for systemic therapy or other adequate infection treatment.
- The interaction of Mepilex Ag with other topical treatments has not been demonstrated.
- Other than saline solution or water, the interaction of cleansing agents in combination with Mepilex Ag has not been demonstrated.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.

Instructions for use

Mepilex Ag should be used by or under the supervision of a qualified health care professional.

Note that local hygiene procedures should be followed prior to and following the dressing change.

1. Cleanse the wound with saline solution or water according to standard clinical practice.
2. Dry the surrounding skin thoroughly.
3. Select appropriate dressing size. For best clinical outcomes, ensure that the product covers the wound bed and overlaps the wound edges.
4. Cut to appropriate size (if needed). Remove the release films and apply the adherent side to the wound. Do not stretch.
5. When necessary, fixate Mepilex Ag with a bandage or other fixation.

Each Mepilex Ag dressing may be left in place for up to 7 days depending on the patient, the condition of the wound and surrounding skin, or as indicated by accepted clinical practice.

Mepilex Ag is intended for short-term use up to 4 weeks. For long-term use, a clinical assessment by a physician is recommended.

Special storage conditions and handling conditions

Mepilex Ag should be stored in dry conditions below 25°C/77°F and protected from direct sunlight. Note that any color variances in Mepilex Ag do not affect product effectiveness or safety.

Disposal should be handled according to local environmental procedures.

Other information

If any serious incident has occurred in relation to the use of Mepilex Ag, you should report to Mölnlycke Health Care and to your local competent authority.

Mepilex and Safetac are registered trademarks of Mölnlycke Health Care AB.

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Mepilex y Safetac son marcas registradas de Mölnlycke Health Care AB.