Mölnlycke®
Melgisorb® Plus
Melgisorb® Plus Cavity

STERILE R




UK CA 0086

CE 2797




Manufacturer / Fabricant
Mölnlycke Health Care AB
Entreprenörsstråket 21, SE-431 53 Mölndal, Sweden

www.molnlycke.com

Melgisorb® Plus

Melgisorb® Plus Cavity

www.molnlycke.com/symbols

Single use
Usage unique

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CH REP

Mölnlycke Health Care AG
Brandstrasse 24, 8952 Schlieren

UK Responsible Person:
Mölnlycke Health Care Ltd.,
Medlock Street, Oldham, Lancashire, OL1 3HS, UK

Toll free number :
USA 1-800-882-4582
Canada 1-800-494-5134


Australian sponsor address:
Mölnlycke Health Care Pty. Ltd.
12 Narabang Way, Belrose, NSW 2085 AUSTRALIA

en Do not use if package is damaged
de Bei beschädigter Verp. verwerfen
fr Ne pas utiliser si l'emballage est endommagé
es No utilizar si el envase está dañado
nl Niet gebruiken indien verpakking beschadigd is
sv Skall inte användas om förpackningen är skadad
it Non utilizzare se la confezione è danneggiata
fi Älä käytä jos pakkaus vaurioitunut
pt Não usar se embalagem danificada
da Brug ikke, hvis pakken er beskadiget
el Μην το χρησιμοποιήσετε αν η συσκευασία είναι κατεστραμμένη
pl Nie stosować, jeśli opakowanie jest uszkodzone
cs Nepoužívejte, je-li obal poškozen
hu Ne használja, ha a csomagolás sérült
no Må ikke brukes hvis forpakningen er skadet
sl Ne uporabite, če je ovoj poškodovan
bg Не използвайте, ако опаковката е увредена
ro Nu utilizați dacă ambalajul e deteriorat
sk Nepoužívajte, ak je obal poškodený
tr Ambalaj hasarlıysa kullanmayın
lt Nenaudoti, jei pakuotė pažeista
lv Nelietot, ja iepakojums ir bojāts
et Kahjustatud pakendi korral mitte kasutada
ru Не использ., если упак. поврежд.
hr Ne koristiti ako je pakiranje oštećeno
sr Ne koristiti ako je omot oštećen
ar لا يجوز استخدام الجهاز إذا كانت العبوة تالفة

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Melgisorb® Plus / Melgisorb® Plus Cavity

Absorbent Alginate Dressing / Absorbent Alginate Cavity Dressing

en

Product description

Melgisorb Plus is a soft, highly absorbent, sterile, single use dressing made of Calcium Sodium Alginate fibres. The alginate is hydrophilic and absorbs wound exudate. When absorbing exudate containing sodium ions, the solid alginate fibres quickly form a moist gel, which creates a moist environment beneficial for wound healing. When in contact with blood, Melgisorb Plus accelerates platelet aggregation, enhancing the haemostatic process. The gel can be easily removed from the fragile wound bed by irrigating with 0.9% saline.

Melgisorb Plus is available as flat and as cavity dressings.

Indications for use

Melgisorb Plus is intended for a wide range of moderately to heavily exuding wounds, both infected and non-infected, such as: pressure sores, venous and arterial ulcers, diabetic ulcers, post-operative wounds, dermal lesions and other external wounds inflicted by trauma.

Precautions

Melgisorb Plus is not indicated for dry wounds, third-degree burns or surgical implantation.

- Melgisorb Plus should not be used on heavily bleeding wounds.
- Do not use Melgisorb Plus on patient and/or user with known hypersensitivity to the ingoing materials/components of the product.
- If infection is suspected, follow local routines given by health care professionals.
- Do not reuse. If reused performance of the product may deteriorate, cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.
- The use of Melgisorb Plus is not recommended for infants less than 12 months of age

Instructions for use

en

Melgisorb Plus can be used by lay persons under supervision of health care professionals.

Application:

1. Cleanse/flush wound with 0.9% saline or appropriate solution.
2. Dry the healthy surrounding skin and if needed, protect with a water-resistant cream/ointment.
3. Apply a dry Melgisorb Plus product to the moist wound bed.

 Shallow wounds: Choose the correct size of flat dressing to be able to cover the entire wound.

 Deep wounds and cavities: Cut the appropriate length of the cavity dressing and pack the wound loosely.
4. Cover with an appropriate secondary dressing. The choice depends on the exudate level. Secure with an adhesive dressing or a bandage as appropriate.
5. Compression therapy may be used in conjunction with Melgisorb Plus.

Dressing change and removal:

The dressing change interval may be several days. Change the dressing before it is fully saturated, or as indicated by clinical practice.

Remove Melgisorb Plus by gently cleansing/flushing with 0.9% saline or appropriate solution. Any non-gelled Melgisorb Plus will moisten in contact with the saline.

Disposal should be handled according to local environmental procedures.

Other information

If any serious incident has occurred in relation to the use of Melgisorb Plus it should be reported to Mölnlycke Health Care and to your local competent authority.

Melgisorb® Plus and Melgisorb® Plus Cavity are registered trademarks of Mölnlycke Health Care AB.

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