Mölnlycke®

# Tortoise® Sling system

www.molnlycke.com/tortoise

## Mölnlycke® Tortoise® Sling system
## Mölnlycke® Tortoise® Sling system bariatric

**Lift compatible turning and positioning system with pressure redistribution**

en

Always follow your facility's safe patient handling policies and procedures when repositioning patients.

The Mölnlycke® Tortoise® Sling system requires two or more caregivers. Use as directed by appropriate clinician and please read all of the instructions before use.

If any serious incident has occurred in relation to the use of Mölnlycke Tortoise Sling system, it should be reported to Mölnlycke Health Care and to your local competent authority.

## Product description

The Mölnlycke Tortoise Sling system is a device that helps to redistribute pressure to a patient's body while also acting as an assistive transfer device for users. The handling of the patient can either be performed by connecting the device to a lift or by manual usage. This device comes in two sizes to accommodate different shapes and sizes of patients.

The device consists of two main parts, one pad and one Mölnlycke® Z-Flo™ Fluidized Positioner (referred to as Z-Flo™ positioner). The pad utilizes an internal, low-pressure air bladder which adapts to the size and shape of a patient, helping to provide pressure redistribution. The outside of the pad is comprised of two durable materials, rip-stop and pack cloth, along with reinforced handles that help users to easily turn and position a patient, with or without a lift, while maintaining proper body mechanics. The Z-Flo™ positioner assists in sacral offloading by maintaining neutral body alignment and a desirable therapeutic angle of the patient. The Mölnlycke Tortoise Sling system is designed to remain under the patient for continuous support during the length of care.

## Product overview

⇨ *See illustration 12.*

**Features and materials**

A & B. Body and tail

*Airbladder (polyurethane), ripstop (polyamide), pack cloth (polyamide), non-skid (polyester)*

The body and the tail includes airbladders which help provide pressure redistribution for patients.

The tail is tucked under the body to improve the envelopment. The tucked tail also hinders patients from sliding. The non-skid material on the tail increases the friction towards the bed and minimizes the device from sliding.

The construction helps to reduce the shear forces during transfer.

C. Black flaps

*Pack cloth (polyamide)*

The black flaps can be tucked under the device to add envelopment. When turning a patient, a black flap is also tucked around the Z-Flo™ positioner to secure and lock it in place.

D. Top piece

*Foam (Polyurethane), Ripstop (polyamide), Pack cloth (polyamide)*

The top piece provides support for patient's head during transfer.

E. White handles

*Straps (polypropylene)*

The white handles along the side of the black flaps are mainly made for manual boosting. The white handles nearest the tail are also used for tucking the tail.

F. Main handles

*Straps (polypropylene)*

The main handles allow users to easily turn and position a patient (with or without a lift). Most handles include several loop settings in order to adjust the length of the strap and the position of the patient. When connecting the handles to a lift, the preferred loops are placed in the hanger bar. In order to adjust the position or remove the device, simply detach the loop from the hanger bar.

G. Labels

*Labels (polyethylene)*

The labels include illustrations and information for correct and safe usage.

H. Z-Flo™ positioner

*Film (polyurethane), microspheres, viscous fluid, beads*

The Z-Flo™ positioner assists in sacral offloading by maintaining neutral body alignment and a desired therapeutic angle of the patient.

## Intended Use / Indications for use / Intended user profile

Intended use: Mölnlycke Tortoise Sling system is designed to provide patient turning and repositioning, which may be used as part of a prophylactic therapy to help prevent pressure ulcers.

Indications for use: Patient turning and repositioning, Pressure redistribution, Lateral transfer, Vertical transfer.

Intended user profile: The product is intended to be used by health care professionals and laypersons under the supervision of health care professional.

## Intended patient population

Mölnlycke Tortoise Sling system can be used on adults and adolescents.

## Contraindications

Do not use Mölnlycke Tortoise Sling system with patients who have unstable thoracic, cervical or lumbar fractures.

## Lifetime

The Mölnlycke Tortoise Sling system may be used with one patient for up to 90 days and should be discarded thereafter.

## Weight limits

Max 272 kg (600 lbs) for Mölnlycke Tortoise Sling system

Max 400 kg (880 lbs) for Mölnlycke Tortoise Sling system bariatric

When connecting the device to a lift, please note that the Safe Working Load (SWL) for each product (lift, hanger bar, and sling) should be greater than the patient weight.

## Imaging

The Mölnlycke Tortoise Sling system (including the Z-Flo™ positioner) is MR Safe and can also be used during CT. The pad can be left under the patient during X-ray, but the Z-Flo™ positioner itself should not be present in the imaging field.

## Warnings

- Always perform a patient assessment (described in the section for Preparations before use) before using the device to secure safe handling.
- Only use a hanger bar with open loop attachments. Never use a hanger bar with clip or other attachments.
- The Safe Working Load (SWL) for each product (lift, hanger bar, and sling) should be greater than the patient weight.
- Never leave patients unattended when lifted.
- Always position the tail of the device under patient's knees before transfer with a lift.
- Always perform a "safety lift" (described in the section for Vertical transfer), to secure that handles are connected accurately to hanger bar, before transferring a patient with a lift.
- Ensure that the patient's knees are positioned higher than their buttocks when the patient is elevated with a lift and/or seated in a chair.
- To avoid puncturing, be careful when handling sharp objects in proximity of the Mölnlycke® Z-Flo™ Fluidized Positioner.
- When handling the Mölnlycke® Z-Flo™ Fluidized Positioner, especially during molding, do not use fingernails and avoid wearing rings with raised or sharp edges.

## Precautions

- The use of a linen is required between the patient's skin and the pad.
- When connecting the device to a lift, a large size hanger bar should be used.
- Pay extra attention to head and feet when transferring a patient.
- Patients who are unable to properly bend their knees or rotate their hips should only use the device if the health care provider considers it to be beneficial to the patient.
- Ensure that the patient's neck is not bent when elevated with a lift.
- The user should be aware of all tubes, IV lines, catheters, and other patient equipment so as not to get stuck or damaged when in use.
- Always ensure that bed and chair brakes are locked when handling the device.
- Keep the handles out of imaging area during X-ray.
- Handles should be tucked away during CT and MRI to avoid entangling.
- Do not X-ray through the Mölnlycke® Z-Flo™ Fluidized Positioner.
- No maintenance required. Not user repairable.
- Replace promptly if punctured or damaged.
- If the Mölnlycke® Z-Flo™ Fluidized Positioner is punctured, place it in a sealed plastic bag to avoid spillage of the fluidized material.

## Other Information

- The device is for single patient use only.
- Product provided non-sterile. Do not sterilize.
- Disposal should be handled according to local environmental procedures.

## Cleaning Instructions

Wipe with hospital approved disinfectant or household antimicrobial cleaner/ disinfectant, including bleach wipes only for the pad, according to the label instructions. Extended use of disinfectants on the Mölnlycke® Z-Flo™ Fluidized Positioner may reduce its lifetime.

## Special storage conditions and handling conditions

Mölnlycke Z-Flo™ Fluidized Positioners should be stored below 40°C (104°F).

Long exposure to high temperatures may deteriorate the product prematurely.

## Compatibility

⇨ ***See illustration 13 and illustration 14.***

- The device is compatible with non-AC-powered patient lifts and AC-powered patient lifts.
- The device is only compatible with hanger bars with open loop attachments. Never use a hanger bar with clip or other attachments.
- The device is recommended to be used with 2-point hanger bar (with open loop attachment), but may be used with 4-point hanger bar (with open loop attachment).

## Instructions for use

The users are responsible for ensuring that the combination between the chosen lift, hanger bar, and Mölnlycke Tortoise Sling system is compatible (see section Compatibility and Weight limits) and that the respective Instructions for Use (IFU) is available and reviewed.

## Preparations before use

a) Perform a patient assessment to determine if the size, type, and shape of the device is appropriate for the patient. The patient's weight, physical disabilities, girth, fat tissue composition, and height are important factors to include in the assessment.

Furthermore, evaluate the number of users needed for potential maneuvers with the patient. At least two caregivers are required for each maneuver.

b) Check the expiry date on the package. Keep in mind that the maximum lifetime of 90 days should not exceed the expiry date.

c) Inspect the device and ensure that it is complete and not damaged. Do not use a damaged device.

d) Read the IFU in full before usage.

e) Keep this IFU for future reference.

f) When the product is being used for the first time, mark the date and patient number/ID on the product to facilitate traceability.

## Product and patient alignment

a) Place the pad on the bed with the spine- and pelvis print facing up and with a draw sheet on top. Make sure to align the patient with the center line and ensure that the patient's sacrum is on top of the pelvis printed on the pad. Always confirm that patient's sacrum is aligned with the lateral sacrum line that goes across the middle of the pad. Also ensure that the same sacrum line aligns where the bed bends, at the hinge joint on the bed.

⇨ ***See illustration 1.***

## Tucking

a) On both sides of the patient, locate the white handles underneath the pad nearest to the tail. Next, take the red handles and tuck the tail under the pad towards the patient's sacrum. The black nonskid material on the red striped tail should be tucked under the pad, against the bed fitted sheet. No red prints or red handles should be visible after engaging the tail.

⇨ ***See illustration 2 and illustration 3.***

b) Tuck the black flaps under the pad itself to create an envelopment effect beneficial for pressure redistribution. By tucking the black flaps under the pad, air is displaced creating greater envelopment.

⇨ ***See illustration 4.***

c) Ensure that all straps are tucked under the black flaps.

d) Adjust the angle of the bed to make it comfortable for the patient.

e) Confirm again that patient's sacrum is aligned with the lateral sacrum line of the pad and the hinge joint on the bed (i.e. where the bed bends).

⇨ ***See illustration 5.***

f) Raise side rails of the bed to reduce risk of falling.

## Turning and Positioning

a) Attach the four middle handles on one side of the pad evenly to the hanger bar. It is recommended to use the outer loop on the handles.

If using a 4-point hanger bar, it is recommended to connect the handles on two diagonally placed attachments.

**WARNING. Only use a hanger bar with open loop attachments. Never use a hanger bar with clip or other attachments.**

⇨ ***See illustration 6.***

b) To have full control when turning the patient, the user on the receiving side of the bed must operate the lift. Make sure that the patient is still on the green part of the pad when gently turning the patient. Placing a pillow (or other padding) between the knees of the patient helps to prevent pressure points and increases comfort.

c) Place the Z-Flo™ positioner between the pad and the bed. Use the lateral sacrum line on the back to locate the sacrum. Mold the Z-Flo™ positioner towards the body from the scapula down towards the buttocks. The Z-Flo™ positioner can be molded further towards the body to increase the degree of turn.

d) Allow the patient to ease back onto the pad. This displaces the air inside the pad and helps to redistribute pressure on the lower hip. Detach the handles from the hanger bar.

e) Tuck the black flap around the Z-Flo™ positioner to secure and lock it in place. Tuck the opposite black flap under the pad until desired envelopment effect is achieved.

⇨ ***See illustration 7.***

f) *If a lift is not available, perform the turn manually. The user on the receiving side grasps two or more of the handles on opposite side and then pulls the pad towards themselves.*

## Vertical transfer (including boosting)

When elevating a patient in the device with a lift, the patient should end up in a recumbent position. This is normally enabled by connecting the outer loops of the four middle handles first, and then adjusting the loops for tail and top handles individually. Regardless of this recommendation it is the caregivers' responsibility to ensure that the chosen loop handles provide a safe position for the patient.

Preparations

a) Preferably place the Z-Flo™ positioner between the patients knees to prevent the hip from internal rotation. This will add extra comfort.

b) Let the patient´s arms rest on the chest or abdomen to increase comfort during elevation with a lift.

c) Attach the four middle handles on each side of the pad evenly to the hanger bar. It is recommended to use the outer loop on the handles.

**PRECAUTION. When connecting the device to a lift, a large size hanger bar should be used.**

**WARNING. The Safe Working Load (SWL) for each product (lift, hanger bar, and sling) should be greater than the patient weight.**

**WARNING. Only use a hanger bar with open loop attachments. Never use a hanger bar with clip or other attachments.**

d) Unfold the tail until it reaches the back of the knee. If needed, assist the patient in bending their knees before attaching the knee straps. Taller patients might need the tail to be fully unfolded and in certain cases supported by an extra caregiver. To improve comfort, make sure to prevent the black non-skid surface from coming in direct contact with skin.

**WARNING. Always position the tail of the device under patient's knees before transfer with a lift.**

e) Proceed by attaching the tail and top handles to the hanger bar. If needed, adjust the length of the handles by using the loops to create tension to the straps.

f) Perform a "safety lift":

Initiate the lift and stop when there is tension in the handles and ensure that all loops are secured.

**WARNING. Always perform a "safety lift", to secure that handles are connected accurately to hanger bar, before transferring a patient with a lift.**

⇨ ***See illustration 8.***

Lifting and transfer

g) Lift the patient carefully. It is recommended that one caregiver holds the patient's feet during the lift for additional security.

**WARNING. Ensure that patient's knees are positioned higher than their buttocks when the patient is elevated with a lift.**

**PRECAUTION. Ensure that patient's neck is not bent when elevated with a lift.**

**PRECAUTION. Pay extra attention to head and feet when transferring a patient.**

⇨ ***See illustration 9.***

**h)** Alternative 1.
For transfer upwards in bed (boosting), simply move the lifted patient upwards and align the sacrum line on the back of the pad with the hinge joint of the bed (i.e. where the bed bends).

Alternative 2.
For transfer to another receiving surface (e.g. cardiac chair, recliner, or another bed), prepare by adjusting the inclination of the seat and the backrest to the same angle as the bed. When transferring the patient to a cardiac chair or recliner, it is important that the receiving surface is tilted a bit backwards, especially the seat. Make sure that both supporting surfaces have locked wheels.

⇨ ***See illustration 9.***

i) Lower the patient towards the chosen receiving surface. Use the sacrum line on the back of pad as a reference for easy alignment of the patient with the surface. Gently place the patient comfortably on the receiving surface.

j) When transferring the patient to a cardiac chair or recliner it is important to provide a safe position. If foot support is available, its use is strongly advised to prevent the patient from sliding.

**WARNING. Ensure that the patient's knees are positioned higher than their buttocks when the patient is seated in a chair.**

k) Smooth out the sheet to ensure no wrinkles are present under the patient. Tuck the tail and the black flaps on both sides of the pad to create envelopment effect. Make sure no straps are stuck to the chair or placed under the patient to avoid pressure points.

⇨ ***See illustration 10.***

l) If applicable, remove the Z-Flo™ positioner from between the patients knees.

m) If applicable, confirm that patient's sacrum is aligned with the lateral sacrum line of the pad and the hinge joint on the receiving surface (i.e. where it bends).

n) For transferring the patient back to the bed repeat the steps above for vertical transfer.

*Manual boosting*

*If a patient lift is not available, boosting can be performed manually.*

*a) Flatten the bed to a horizontal position. Locate the white boosting handles on both sides under the black flaps.*

*b) With arms at hip distance and weight evenly distributed, use appropriate ergonomics to boost the patient up in the bed.*

*c) Confirm that patient's sacrum is aligned with the lateral sacrum line of the pad and the hinge joint on the bed (i.e. where the bed bends).*

*d) Tuck the tail and the black flaps. Set the bed in a comfortable position.*

## Lateral transfer

*If there is no access to a patient lift, the patient may be transferred manually. It is recommended to include four to five caregivers for this maneuver.*

⇨ ***See illustration 11.***

*a) Place both surfaces in close proximity to each other. Make sure the bed brakes are locked on both beds. Ensure the surface the patient is being transferred to is at same height or slightly lower than the surface they are being transferred from.*

*b) Untuck the tail for an easier transfer.*

*c) If desired, position a lateral sliding aid between the original surface and the end point surface.*

*d) Position members of the patient transfer team on both sides of the patient to prevent the beds from moving apart. If needed, additional caregivers can hold the patient's feet and head/neck to avoid drag.*

*e) Coordinate timing with transfer team and then gently pull the handles to transfer the patient slow and steady.*

## Troubleshooting

| Problems | Solutions |
|---|---|
| It is hard to reach the hanger bar with the tail and/or top strap. | - Bend the patient's knees to decrease the distance between tail strap and hanger bar.<br>- Adjust the angles in bed. Raise the backrest and elevate the knees.<br>- Chose the bariatric size of the device that has wider girth. |
| Patients with fragile skin experience discomfort in the back of the knees when elevated with a lift. | Place pillows or other padding under the knees when lifting. |
| It is not possible to fold the tail behind the patient's knees. This can typically happen for very tall patients. | Untuck the tail completely. If the edge of the tail and/or non-skid surface is causing discomfort for the patient place pillows or other padding in between. |
| The patient complains that it is uncomfortable to be elevated using a lift. | - Adjust the loop of the handles. Aim to enable a recumbent position and avoid that the patient's neck is bent.<br>- Place a Z-Flo™ positioner between the patient's knees to prevent the hip from internal rotation during elevation with a lift.<br>- Let the patient´s arms rest on the chest or abdomen to increase comfort during elevation with a lift. |

Additional tips and information can be found using the QR code.

www.molnlycke.com/symbols

Caution, see instructions for use

Consult instructions for use

Medical Device

MR Safe

Upper temperature limit

Patient number

Write

First day in use

Wipe off

Do not wash

Do not bleach

Do not tumble dry

Do not iron

Do not dry clean

Body weight

Tortoise Sling: Max 272 kg / 600 lbs
Tortoise Sling bariatric: Max 400 kg / 880 lbs

Australian sponsor address:
Mölnlycke Health Care Pty. Ltd.
12 Narabang Way, Belrose, NSW 2085, AUSTRALIA

Patent number :
US 8.171.585 B2 US 9.504.621 B2
US 9.120.666 B2 US 9.833.371 B2

Manufacturer / Fabricant
Mölnlycke Health Care AB
Entreprenörsstråket 21, SE-431 53 Mölndal, Sweden
www.molnlycke.com

The Mölnlycke trademark, name and logo are registered globally to one or more of the Mölnlycke Health Care group of companies

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