Tubigrip®

High-quality[1] multi-purpose elasticated tubular bandage

Firm support for sprains, strains and swelling

Discover the benefits of Tubigrip

- Easy to apply – no need for pins or tape[1]
- Provides continuous, uniform support
- Durable – can be washed and reused
- Cut to length for a custom fit
- Multiple sizes to suit most applications
- Provides light compression[2–4]
- Can be used in the management of oedema[2,4]



Find the correct size

Tubigrip size	Minimum limb circumference (cm)	Minimum limb circumference (inch)	Maximum limb circumference (cm)	Maximum limb circumference (inch)	Sub-bandage pressure (mmHg)
A	10	3.9	14.5	5.7	20–30
B	14.5	5.7	29	11.4	10–20
C	22	8.7	35	13.8	
D	26	10.2	40	15.7	
E	29	11.4	46	18.1	
F	35	13.8	53	20.9	
G	40	15.7	62	24.4	5–15
J	54	21.3	82	32.3	
K	70	27.6	113	44.5	
L	106	41.7	149	58.7	
M	119	46.9	167	65.7	

The sub-bandage pressures are obtained in double layer.

1. Mönlycke Health Care. Data on file. **2.** Gleeson AP, Stuart MJ, Wilson B, Phillips B.Ultrasound assessment and conservative management of inversion injuries of the ankle in children: Plaster of Paris versus Tubigrip. Journal of Bone and Joint Surgery - Series B.1996;78(3):484–487. **3.** Naeem M, Rahimnajjad MK, Rahimnajjad NA, Idrees Z, Shah GA, Abbas G. Assessment of functional treatment versus plaster of Paris in the treatment of grade 1 and 2 lateral ankle sprains. Journal of Orthopaedics and Traumatology. 2014. **4.** Tivendale L, Wynne R, Law SS, Kenny L, Dean G, Tatoulis J. Effect of external compression bandage on forearm swelling, discomfort, and pain following the removal of a radial artery conduit for CABGS. Chest.2011;140(4).

Find out more at molnlycke.com

Mölnlycke Health Care Limited, Unity House, Medlock Street, Oldham, Lancashire, OL1 3 HS, United Kingdom. Tel: +44 (0) 161 621 3900
The Mölnlycke, Mepilex, trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2025 Mölnlycke Health Care Limited. All rights reserved. UKWC1358

Mölnlycke®